


Obstruction Removal Kit for ST Nebulizer Nebulizers

On occasion, the ST nebulizers may become clogged. This kit is specifically designed for PFA-ST and PolyPro-ST nebulizers.

- Includes a guide tool and 10 pieces of fused silica.
- 0-3 ml/min sample flow rate with peristaltic pumping.
- Easy to use:
 - tighten guide onto back of ST nebulizer
 - slide fused silica through
- Silica refills also available.



Obstruction Removal Kit shown with PFA-ST Nebulizer (not included)



Ordering Information

ES-2046-0001	Fused Silica Obstruction Removal Kit for ST type Nebulizers. Includes one guide, 10 pieces of fused silica and instructions.
ES-2046-0002	Fused Silica Refill for Obstruction Removal Kit For use with ST type Nebulizers only. Package of 10.